

V-TRAK

pico

Dynamic support for kids



Inspired by the therapist - The integral lateral support was inspired from experiencing clinic after clinic. With the therapist using their hands to show where they wanted support. This support was designed into our backrests to replicate hands supporting the child



It's what's inside that counts - Throughout the V-Trak range of backs all the adjustment is in the backrest.



Because Kids will be kids - Growth movement and durability. The pico range is designed for kids. The adjustability in the backrest allows for changes to be made easily for weight gain and growth.



Backrest Extension - BKXT-S



Can be inverted to provide thoracic or pelvic support



Accessories - Universal Harness carrier



VP-DC1 - Deep contour



VP-P1 - Profile Shallow contour

Technical Data

| | | Height (cm) | Outside Width (cm) | Inside / User Width (cm) | Contour Depth (cm) | Weight (kg) |
|---------------|---|-------------|--------------------|--------------------------|--------------------|-------------|
| VP-P1 | V-Trak Paediatric profile backrest | 24 | 30 > 33 | 23 > 26 | 5 > 7 | 0.8 kg |
| VP-DC1 | V-Trak Paediatric Deep contour backrest | 24 | 32 > 38 | 26 > 32 | 7 > 10 | 0.9 kg |
| BKXT-S | Backrest extension | 10 | 18 | N/A | N/A | 0.17 kg |

All measurements are an approximate as the foam will compress and should be used as a guide